SELECIONE O IDIOMA

SELECT THE LANGUAGE

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
Intraoral Aesthetic Archwire

Intraoral Aesthetic arches in NiTi Superelastic are components of the orthodontic appliance, designed to transmit light clinical forces generated by the interaction between the arch and brackets, promoting the movement of dental elements.

They can easily adapt to the device even in cases of severe misalignment thanks to its ability to deform by light force action and return to its original shape, this feature allows the Intraoral Arch in Super Elastic NiTi to generate smooth and constant forces throughout the movement without losing its original shape.

The aesthetic intraoral arches have a similar color to the dental enamel, which presents an excellent aesthetic result and leaves the device more discreet.

Use Indication:

The basic objective of an intraoral arch is to achieve optimal occlusion, with well aligned arches.


The most common clinical applications are:

- a) Angle's Classification:
 - Class I;
 - Class II;
 - Class III.
- b) Problems on teeth and bones:
 - Prognathism;
 - Retrognathism;
 - Deep bite (overbite);
 - Openbite;

- Crowding (lack of spaces);
- Anterior crossbite;
- Posterior crossbite;
- Diastema.

Identification of archs:

The arches must be chosen according to the characteristics of the clinical case, they are identified by the markings of medium II or large size III.


Deflection curves:

During orthodontic treatment, the arches will undergo temporary deformation due to dental misalignment. These deformations are classified according to the direction of the forces involved:

- Lingual-Vestibular;
- Occlusal-Gingival.

The efforts exerted during the treatment should be estimated considering the measurements made in the deflection plateau, which represents the behavior of the wire during clinical use (on the chart the line indicated by "a").

The typical behavior of a Type 2 arch is shown in the graph below, according to the test provided in EN ISO 15841, typical and obtained values are shown in sequence.

Legend: X deflection, in mm; Y force, in Newtons; the results are measured in the deflection plateau (bottom curve).

1 Newton = 102 gf

Deflection points	Ø .012"	Ø .014"	Ø .016"	Ø .018"
0,5 mm	0,27 N	0,59 N	1,20 N	1,90 N
1 mm	0,57 N	1,11 N	1,69 N	2,44 N
2 mm	0,59 N	1,19 N	1,76 N	2,66 N
3 mm	0,69 N	1,21 N	1,77 N	2,90 N

Direction of force application - Lingual-Vestibular		
Deflection points	.017" x .025"	.019" x .025"
0,5 mm	4,80 N	5,57 N
1 mm	5,58 N	6,13 N
2 mm	6,83 N	7,70 N
3 mm	8,60 N	9,83 N

Direction of force application - Occlusal-Gingival		
Deflection points	.017" x .025"	.019" x .025"
0,5 mm	3,48 N	4,42 N
1 mm	4,14 N	4,80 N

2 mm	4,48 N	5,59 N
3 mm	5,10 N	6,62 N

Type 2 wire according to EN ISO 15841, because it has non-linear deflection.

Contraindications:

Morelli supplies the products only to qualified professionals. It is the responsibility of the orthodontist to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene;
- Patients unable to cooperate with treatment;
- Patients with propitious oral environment to dental enamel demineralization;
- Allergy to any of the materials of the orthodontic appliance;
- Diseases or other pre-existing conditions that may hinder orthodontic treatment;
- Root resorption;
- Existing bone resorption;
- Existing Decalcification in dental enamel;
- Use in a Magnetic Resonance environment.

Side effects:

It is the responsibility of the qualified professional to identify any possibility of risk of harm and / or contraindication that may arise during treatment, communicate to the patient the possibility of any unwanted side effect.

During treatment, unwanted side effects may include: damage to oral mucosa and tongue, difficulty in speaking or chewing, tooth discoloration, decalcification, root resorption, periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, discomfort, pain, discomfort and sensibility.

Composition/ Information of components:

Ni = 55,4% a 57% / Ti = Bal.

Composition of Coating:

Low cure thermofix urethane copolymer resin (Polyurethane).

Expiry date:

10 years

Warnings:

Thermal treatments are not recommended.

Avoid handling the arch with pliers or tweezers, and install accessories such as stops and hooks carefully.

The durability of the aesthetic layer varies with the oral environment: Consider replacing them within 30 days or earlier to maintain the desired aesthetic effect.

Heating the arc above 200 ° C can induce fracture and deteriorate the aesthetic layer. Overstressed and sharp bends can fracture the arch.

Frequent exposure to substances with fluoride may weaken the arch.

For distal cutting in the buccal environment, make sure that the pliers retain the tip of the Arch.

This product contains Nickel. A small percentage of the population is allergic to these metal. In case of allergic reactions, advise the patient to seek medical attention.

Frequent exposure to substances such as fluoride and alcohol can cause the Arch to weaken and decrease the resistance of the aesthetic layer.

This product must not be used in a Magnetic Resonance environment as it may cause interference with the images during the examination. It is up to orthodontic practitioner to warn the patient about this contraindication and arrange for removal of the device if necessary.

Non reusable product, as its reuse may cause cross infection and loss of mechanical properties due to natural wear.

The product is supplied in unsterile condition, however has proven biocompatibility, not being necessary cleaning, disinfection or sterilization measures, since the package has not been violated. If the product comes into contact with surfaces or substances other than its intended purpose prior to use, it is recommended to discard the product. The product should not be reused or reprocessed.

Product for exclusive application by a qualified professional, the use by person without the necessary technical knowledge may cause unwanted tooth movement, bone loss, loss of the dental element, fenestration, gingival recession or root dehiscence.


It is up to the qualified professional to guide the patient about correct oral hygiene, to avoid the appearance of plaque and tartar, or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis.

Considering that the product is used in contact with mucous membranes and body fluids, it is recommended to the professional to apply the sanitary norms applicable for the disposal of medical products according to current sanitary regulations.

ATTENTION:

Morelli is not responsible for the results obtained by applying the product without due precautions or non-observance of the warnings.

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