SELECIONE O IDIOMA

SELECT THE LANGUAGE

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
Orthodontic Wires (Type I)*

* Type I wire according to EN ISO 15841, by having linear elastic behavior

The wires are fixed in the slots of brackets, tubes or other orthodontic accessories and through the unevenness of the teeth correct the positioning of the teeth based on an ideal pre-contoured shape of the dental arch to be followed. The wires can also be applied for fixing and anchoring orthodontic accessories.

Use Indication:

Orthodontic assistant indicated for the treatment of class I, II or III malocclusions, correcting teeth and bone problems such as:

- Prognathism;
- Retrognathism;
- Deep bite;
- Open bite;
- Crowding;
- Anterior crossbite;
- Posterior crossbite.

Contraindications:

Morelli supplies the products to qualified professionals only. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with oral hygiene deficiency;
- Patients unable to cooperate with treatment;
- Allergy to any of the orthodontic appliance materials;
- Diseases or other pre-existing conditions that may hinder orthodontic treatment;
- Pronounced root resorption;

- Existing bone resorption;
- Use in magnetic resonance environment.

Side Effects:

It is the orthodontist's sole responsibility to identify any possible risk of harm and / or contraindication that may arise during treatment, to report to the patient the possibility of any unwanted side effects.

During treatment, unwanted side effects may include: damage to the oral mucosa and tongue, difficulty speaking or chewing, tooth discoloration, decalcification, root resorption, periodontal complications, allergic reactions, difficulty maintaining oral hygiene, pain, discomfort and sensitivity.

Composition/ information about components:

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10 years

Warnings:

Thermal treatments are not recommended.

Frequent exposure to substances with fluoride may weaken the arch.

For distal cutting in the buccal environment, make sure that the pliers retain the tip of the Arch.

This product must not be used in a Magnetic Resonance environment as it may cause interference with the images during the examination. It is up to orthodontic practitioner to warn the patient about this contraindication and arrange for removal of the device if necessary.

This product contains nickel and chrome. A small percentage of the population is allergic to these metals. If allergic reactions occur, advise the patient to seek medical advice.

Non reusable product, as its reuse may cause cross infection and loss of mechanical properties due to natural wear.

The product is supplied in unsterile condition, however has proven biocompatibility, not being necessary cleaning, disinfection or sterilization measures, since the package has not been violated. If the product comes into contact with surfaces or substances other than its intended purpose prior to use, it is recommended to discard the product. The product should not be reused or reprocessed.

It is up to the qualified professional to advise their patients of the risk of aspiration and swallowing in the event of breakage / detachment of parts, should not underestimate the possible immediate or late complications resulting from this type of accident and need to adequately inform the patient and / or their legal guardians about the risks inherent in the event and the actions to be taken.

Do not use the product if its packaging is damaged or opened in order to minimize possible dangers related to the transmission of microorganisms or possible falsification of products.

Product for exclusive application by a qualified professional, the use by person without the necessary technical knowledge may cause unwanted tooth movement, bone loss, loss of the dental element, fenestration, gingival recession or root dehiscence.


It is up to the qualified professional to guide the patient about correct oral hygiene, to avoid the appearance of plaque and tartar, or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis.

Considering that the product is used in contact with mucous membranes and body fluids, it is recommended to the professional to apply the sanitary norms applicable for the disposal of medical products according to current sanitary regulations.

ATTENTION:

Morelli is not responsible for the results obtained by applying the product without due precautions or non-observance of the warnings.

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