SELECIONE O IDIOMA

SELECT THE LANGUAGE

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## Intraoral Archwires Copper NiTi

Intraoral arches in Copper NiTi are components of the orthodontic appliance, designed to transmit light clinical forces generated by the interaction between the arch and brackets, promoting the movement of dental elements.

They can easily adapt to the device even in cases of severe misalignment due to its ability to be deformed by slight force action and return to its original shape, this feature allows the intraoral Copper NiTi Arches generate smooth and constant forces throughout the orthodontic movement without losing its original form.

### Use Indication:

When cooled the arches become malleable, allowing the use in severe misalignments, the elevation to the activation temperature causes the arc to return to its original shape, providing:

- Constant unloading plateau, guaranteeing the generation of forces that are biologically favorable to dental movements for longer periods of time.
- More comfort during treatment;
- Greater resistance to permanent deformation.

Such characteristics allow a more consistent treatment planning, due to the greater stability in relation to force versus deformation.

The basic objective of an intraoral arch is to achieve optimal occlusion, with well aligned arches.

The most common clinical applications are:

#### a) Angle's Classification:


- Class I;
- Class II;
- Class III.

#### b) Problems on teeth and bones:

- Prognathism;
- Retrognathism;
- Deep bite (overbite);
- Openbite;
- Crowding (lack of spaces);
- Anterior crossbite;
- Posterior crossbite;
- Diastema.

### Identification of archs:

The arches must be chosen according to the characteristics of the clinical case, they are identified by the markings of lower arch II or upper arch III.

### Mechanical properties

#### Final Austenitization Temperature (Af)


Intraoral arches in Copper NiTi have a final activation temperature of 35°C.

### Deflection curves:

During orthodontic treatment, the arches will undergo temporary deformation due to dental misalignment.

The efforts exerted during the treatment should be estimated considering the measurements made in the deflection plateau, which represents the behavior of the wire during clinical use (on the chart the line indicated by "a").

The typical behavior of a Type 2 arch is shown in the graph below, according to the test provided in EN ISO 15841, typical and obtained values are shown in sequence.

Legend: X deflection, in mm; Y force, in Newtons; the results are measured in the deflection plateau (bottom curve).

1 Newton = 102 gf

Profile	Dimensions	Deflection Force Rate FA N/mm				
		Minimal Diagonal	2º Ordem - (occlusion-gingival)			
			0,5mm	0,5mm	0,5mm	0,5mm
○	∅ .014" - 0,35mm	-	0,59	0,61	0,53	0,63
	∅ .016" - 0,40mm	-	1,0	1,05	0,98	1,031
	∅ .018" - 0,45mm	-	0,98	1,11	1,29	1,69
□	.014" x .025" - 0,35x0,63mm	0,655mm	1,08	1,16	1,26	1,65
	.016" x .022" - 0,40 x 0,55mm	0,617mm	1,67	1,79	1,88	2,11
	.017" x .025" - 0,43 x 0,63mm	0,693mm	1,60	1,78	1,97	2,76
	.018" x .025" - 0,45 x 0,63mm	0,709mm	1,81	2,00	2,21	2,92
	.019" x .025" - 0,48 x 0,63mm	0,724mm	1,63	1,89	2,41	3,64

**Contraindications:**

Morelli supplies the products only to qualified professionals. It is the responsibility of the orthodontist to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene;
- Patients unable to cooperate with treatment;
- Patients with propitious oral environment to dental enamel demineralization;
- Allergy to any of the materials of the orthodontic appliance;
- Diseases or other pre-existing conditions that may hinder orthodontic treatment;
- Root resorption;
- Existing bone resorption;

- Existing Decalcification in dental enamel;
- Use in a Magnetic Resonance environment.

**Side effects:**

It is the responsibility of the qualified professional to identify any possibility of risk of harm and / or contraindication that may arise during treatment, communicate to the patient the possibility of any unwanted side effect.

During treatment, unwanted side effects may include: damage to oral mucosa and tongue, difficulty in speaking or chewing, tooth discoloration, decalcification, root resorption, periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, discomfort, pain, discomfort and sensibility.

**Composition/ Information of components:**

Ni = 49% / Cu= 6% / Ti = Bal.

**Expiry date:**

10 years

**Warnings:**

Thermal treatments are not recommended.

Heating the arch above 200 ° C can induce fracture.

Excessive stresses and sharp folds can fracture the arch.

Frequent exposure to substances with fluoride may cause the weakening of the Arch.

For distal cutting in the buccal environment, make sure that the pliers retain the tip of the Bow.

This product must not be used in a Magnetic Resonance environment as it may cause interference with the images during the examination. It is up to orthodontic practitioner to warn the patient about this contraindication and arrange for removal of the device if necessary.

This product contains Nickel. A small percentage of the population is allergic to these metal. In case of allergic reactions, advise the patient to seek medical attention.

Non reusable product, as its reuse may cause cross infection and loss of mechanical properties due to natural wear.

The product is supplied in unsterile condition, however has proven biocompatibility, not being necessary cleaning, disinfection or sterilization measures, since the package has not been violated. If the product comes into contact with surfaces or substances other than its intended purpose prior to use, it is recommended to discard the product. The product should not be reused or reprocessed.

It is up to the qualified professional to guide their patients from the risk of aspiration and deglutition in the hypothesis of breakage / detachment of the pieces, should not be underestimated the possible immediate or late complications resulting from this type of accident and must adequately inform the patient and / or their legal guardians about the risks inherent in the event and what actions to take.

Do not use the product if its packaging has been damaged or opened, in order to minimize possible dangers related to the transmission of microorganisms or possible falsification of products.

Product for exclusive application by a qualified professional, the use by person without the necessary technical knowledge may cause unwanted tooth movement, bone loss, loss of the dental element, fenestration, gingival recession or root dehiscence.

It is up to the qualified professional to guide the patient about correct oral hygiene, to avoid the appearance of plaque and tartar, or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis.

Considering that the product is used in contact with mucous membranes and body fluids, it is recommended to the professional to apply the sanitary norms applicable for the disposal of medical products according to current sanitary regulations.

**ATTENTION:**

Morelli is not responsible for the results obtained by applying the product without due precautions or non-observance of the warnings.

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