SELECIONE O IDIOMA

SELECT THE LANGUAGE

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## Extrarradicular Orthodontic Steel Mini Screw

Mini Screws are single-use implantable medical and dental devices, temporarily fixed to the maxilla or mandible bone in extra-alveolar areas such as the Mandibular Shelf and Infracrestal Crest to create stable anchor points that are unable to move by applying the clinical forces necessary for orthodontic treatment.

This anchorage, called skeletal anchorage, allows the application of continuous forces and replaces traditional anchor systems or those that require patient collaboration, as well as enabling more complex orthodontic movements with precision and in less time.

**Use indication:**

Description of the most common clinical application situations:

- Alveolar distalization of the entire arch in class II patients;
- Anchorage for block retraction in cases of extraction;
- Bi protrusion treatment;
- Canine and premolar distalization;
- Mesialization of molars;
- Mesialization with protrusion of all teeth in patients with maxillary deficiency or mandibular retrusion;
- Intrusion of posterior teeth;
- Correction of occlusal plane asymmetry;
- Correction of midline deviations;

- Traction of canines included;
- Upper arch retraction for preparation of class III orthognathic surgery;
- Early treatment of class III;
- Brodie cross bite correction.

Mini-screws are preferably recommended for patients under the following conditions:

- Non-cooperative;
- With a reduced number of dental elements;
- In need of maximum anchorage;
- In need of complex movements if traditional anchoring methods are applied.

#### Installation and Removal:

Proper installation tools designed specifically for this application should be used for installation and removal.

The **"Toolkit"** \* (ref.37.20.001) consists of Hand Key with Interchangeable Tips that attach to the hexagonal head portion of the mini-screw, allowing for proper installation or removal. The Interchangeable Tips, besides allowing mounting on the Hand Key, have standard fitting for surgical contra-angle.

\* All parts of the Toolkit are autoclavable..

#### Selection of insertion location:

Extraradicular Mini Screws are preferably installed in regions of high bone density, parallel to the midline of the teeth avoiding root contact.

In the mandible, it is recommended to have the Mini Screws in the region of greater bone margin, present between the first and third lower molars in the region called buccal shelf.

In the maxilla it is recommended to install the bony protuberance along the curvature between the alveolar process and the zygomatic bone, in the region called the infrazigomatic crest.

#### Contraindications:

The patient applying for the mini-screw installation should go through an anamnesis to verify possible contraindications to the use of this anchoring method. Patients with systemic diseases such as type I diabetes, anemia, weakened immune system, local and systemic bone disorders, patients undergoing radiotherapy, poor oral hygiene and pregnant patients should be carefully evaluated due to the stress involved and the possibility of gingivitis gravidarum. In the presence of any of these conditions, a medical evaluation accompanied by the orthodontist is required to make sure that all implications are understood.

#### Side effects:

It is the orthodontist's sole responsibility to identify any potential risk of injury and / or contraindication that may arise during treatment, to report to the patient the possibility of any unwanted side effects.

During treatment, unwanted side effects may include: damage to the oral mucosa and tongue, difficulty speaking or chewing, root resorption, periodontal complications, allergic reactions, difficulties maintaining oral hygiene, pain, discomfort and tenderness. Bone malformation or infection that leads to loosening and loss of stability of the Mini Screw.

This surgical procedure can not only cause the side effects and complications mentioned above, but also problems such as nerve damage, infections, pains that are not necessarily caused by the implant.

#### Composition / information about the components:

Stainless steel

(Cr = 15% - 17,5% / Ni = 3,0% - 5,0% / Cu = 3% - 5% / Mn = Máx 1,5% / Si = Máx 0,7% / Nb = 0,15% - 0,45%)

#### Expiry date:

3 years

#### Warnings:

This product contains nickel and chrome. A small percentage of the population is allergic to these metals. If allergic reactions occur, advise the patient to seek medical advice.

Only open the sterilized envelope at the time of installation of the Mini Screw to ensure sterilization until the time of use.

Do not use the Mini Screw if it has expired or the envelope has been tampered with.

Do not subject the Mini Screw to forces greater than 450gf.

Do not exceed the maximum torque of 55N.cm during installation. Literature reports that fractures may occur due to flexion efforts during application. It is up to the practitioner to decide whether to remove the implant fragment or install another Mini Screw.

In case of immediate loading on this device, the application should be done perpendicular to the long axis of the Mini Screw.

In case of high bone density, use drill up to Ø 1.5 mm for pre-drilling.

Keep the minimum distance of 1mm between the dental roots and the Mini Screw.

Observe the direction of insertion of the Mini Screw to avoid unscrewing while applying mechanics with orthodontic wires.

This product should not be used in an MRI environment as it may cause interference with images during the exam. It is up to the orthodontic professional to warn the patient about this contraindication and to arrange the removal of the device if necessary.

Product for exclusive application by a qualified professional, the use by person without the necessary technical knowledge may cause unwanted tooth movement, bone loss, loss of the dental element or even irreversible damages.

It is up to the qualified professional to guide the patient about correct oral hygiene, to avoid the appearance of plaque or diseases such as gingivitis, periodontitis or even endocarditis.

Considering that the product is used in contact with mucous membranes and body fluids, it is recommended to the professional to apply the sanitary norms applicable for the disposal of medical products according to current sanitary regulations.

**ATTENTION:**

Morelli is not responsible for the results obtained by applying the product without due precautions or non-observance of the warnings.

Do not resterilize the Mini Screw.

Not reusable product. Sterility, cleanliness and mechanical strength of the device are not guaranteed in case of reuse.

Do not use the product if its packaging is violated to minimize possible dangers related to the transmission of microorganisms or possible falsification of products.

Product sterilized by ethylene oxide. Because it is used invasively in the human body can present risks associated with its use, so its sterility must be guaranteed.

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