SELECIONE O IDIOMA

SELECT THE LANGUAGE

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Rhodium Superelastic Intraoral Arch – NiTi

Rhodium Superelastic Intraoral Arches are components of the orthodontic appliance, designed to transmit light clinical forces generated by the interaction between the arch and the brackets, promoting the movement of the dental elements.

They can easily adapt to the device even in cases of severe misalignment thanks to its ability to deform by light force action and return to its original shape, this feature allows Rhodium Superelastic Intraoral Arcs to generate smooth and constant forces throughout orthodontic movement without losing its original form.

Rhodium applied on the surface of the arch, promotes and enhances the aesthetic effect ensuring a discreet touch for the treatment.

Use Indication:

The basic goal of an intraoral arch is to achieve optimal occlusion with well-aligned arches. The most common clinical applications are:


- a) Angle's Classification:
 - Class I;
 - Class II;
 - Class III.

- b) Problems in teeth and bones:
 - Prognathism;
 - Retrognathism;

- Deep Overbite;
- Open bite;
- Crowding;
- Anterior crossbite;
- Posterior crossbite;
- Diastema.

Identification of arches:

The Rhodium Superelastic Intraoral Arches – NiTi should be chosen according to the characteristics of the clinical case, they are identified by the lower II or upper III arch markings.


Deflection curves:

During orthodontic treatment, the arches will undergo temporary deformation due to dental misalignment. These deformations are classified according to the direction of the forces involved:

- Lingual-Vestibular;
- Occlusal-Gingival.

The efforts exerted during the treatment should be estimated considering the measurements made in the deflection plateau, which represents the behavior of the wire during clinical use (on the chart the line indicated by "a").

The typical behavior of a Type 2 arc is shown in the chart below, according to the test provided for in EN ISO 15841, in sequence, typical and obtained values are shown.

Legend: X deflection, in mm; Y force, in Newtons; the results are measured in the deflection plateau (bottom curve).

1 Newton = 102 gf

Deflexion Points	Ø .012"	Ø .014"	Ø .016"	Ø .018"
0,5 mm	0,33	0,89	1,56	2,16
1 mm	0,49	1,11	1,72	2,30
2 mm	0,46	1,09	1,69	2,38
3 mm	0,42	1,01	1,73	2,38

Direction of force application – Lingual-Vestibular		
Deflexion Points	.017" x .025"	.019" x .025"
0,5 mm	5,25	6,15
1 mm	6,02	7,10
2 mm	7,10	8,52
3 mm	7,38	9,00

Direction of force application – Occlusal-Gingival		
Deflexion Points	.017" x .025"	.019" x .025"
0,5 mm	4,20	5,41

1 mm	4,23	5,73
2 mm	4,49	6,40
3 mm	4,66	6,64

Type 2 wire according to EN ISO 15841, because it has non-linear deflection.

Contraindications:

Morelli Ortodontia supplies the products only to qualified professionals. It is the responsibility of the orthodontist to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral care;
- Patients incapable of cooperating with the treatment;
- Allergy to any of the orthodontic device materials;
- Diseases or other pre-existent conditions which may impair orthodontic treatment;
- Existing bone resorption;
- Existing Decalcification in dental enamel.

Side effects:

It is the responsibility of the qualified professional to identify any possibility of risk of harm and / or contraindication that may arise during treatment, communicate to the patient the possibility of any unwanted side effect.

During treatment, unwanted side effects may include: tooth discoloration, decalcification, root resorption, periodontal complications, allergic reactions, difficulties in maintaining oral hygiene, discomfort and pain.

Composition/ Information of components:

Ni = 55,4% a 57% / Ti = Bal.

Expiring date:

10 anos

Warnings:

Thermal treatments are not recommended.

Avoid damage to the aesthetic layer, replace the intraoral arch within 30 days.

Avoid handling the arch with pliers or tweezers, and install accessories such as stops and hooks carefully.

Heating the arc above 200 ° C can induce fracture and deteriorate the aesthetic layer.
Overstressed and sharp bends can fracture the arch.

Frequent exposure to substances with fluoride may weaken the arch.

For distal cutting in the buccal environment, make sure that the pliers retain the tip of the Arch.

This product hasn't been evaluated for compatibility and safety in Magnetic Resonance environment, the patient should be advised of this condition.

This product contains Nickel. A small percentage of the population is allergic to these metal. In case of allergic reactions, advise the patient to seek medical attention.

Non reusable product, as its reuse may cause cross infection and loss of mechanical properties due to natural wear.

The product is supplied in unsterile condition, however has proven biocompatibility, not being necessary cleaning, disinfection or sterilization measures, since the package has not been violated. If the product comes into contact with surfaces or substances other than its intended purpose prior to use, it is recommended to discard the product. The product should not be reused or reprocessed.

Product for exclusive application by a qualified professional, the use by person without the necessary technical knowledge may cause unwanted tooth movement, bone loss, loss of the dental element or even irreversible damages.

It is up to the qualified professional to guide the patient about correct oral hygiene, to avoid the appearance of plaque or diseases such as gingivitis, periodontitis or even endocarditis.

Considering that the product is used in contact with mucous membranes and body fluids, it is recommended to the professional to apply the sanitary norms applicable for the disposal of medical products according to current sanitary regulations.

ATTENTION:

Morelli Ortodontia is not responsible for the results obtained by applying the product without due precautions or non-observance of the warnings.

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